

CATEGORY SPOTLIGHT



VIDEO

See NoCord VPS used in a clinical case involving PFM crown preparation on maxillary centrals.



Taking Impressions with Confidence

NoCord VPS is a single-step, self-retracting impressioning system designed for ease of use and greater predictability

Show up to a job interview, and the first impression can often be make-or-break. Nice attire, fresh breath, and a well-put-together appearance can go a long way in relaying an image of confidence and success. When making a first impression in the literal context of dentistry, the confidence should be within the impression material itself, as successful restorative treatment depends on the material's accuracy, detail, and clarity.

Centrix's NoCord VPS is a one-step, self-retracting impressioning system designed specifically for the wash and tray technique that consists of the Centrix NoCord Wash and NoCord MegaBody Tray Materials.

Without a separate retraction procedure, NoCord Wash and MegaBody Tray Materials have been specially formulated to work together as a system that retracts, controls bleeding, and takes a clear, final impression, while eliminating the time and cost of using retraction cords or pastes and applying an astringent.

"I had been looking for an impression material that did not harm sensitive gum tissues. I work with patients and use my laser to get their gums exactly where I want, so they have their best smile moving forward. NoCord fits that bill," said speaker and lecturer Dr. Nancy Summer Lerch.

NoCord Wash is the first wash material to incorporate a hemostatic agent, a 15% alum (ammonium aluminum sulfate), that effectively closes cut capillaries, controlling bleeding and fluid flows. When paired with the "stiffer" NoCord



MegaBody Tray Material, the wash is driven into the sulcus to displace gingiva and completely capture prep margins. No retraction cord means no gingival trauma for the patient and a less time-consuming process for the clinician.

"NoCord VPS also supplies hemostasis and removes the need for extra materials at the time of impressions. On further use, I was happy to discover the ease of application and the ideal fit of the final restorations from my esthetic dental lab," Dr. Lerch added.

TOP 5 REASONS WHY DR. LERCH USES NOCORD VPS:

- 1 Easy to use
- 2 Atraumatic to the tissues
- 3 Provides hemostasis
- 4 Very visible, easy-to-read margins
- 5 Consistency—works every time



FOR FREE INFORMATION:

800.921.4806
ext. 181

www.dps.li/a/5R9-240

Circle 240
on the card